


Single-sterile/BNS Alternative Kitter Affected Lot FAQ

How do I know if my customer is affected?

- Only the affected lot numbers are impacted. Therefore, if your customer is utilizing a material listed below and it **does not** have the affected lot number format they **may continue to utilize** the product.
- If the customer is utilizing a material listed in Table 1 on the following page with the associated affected lot number format, they must discontinue use and segregate from current inventory. Further instructions will be provided as part of the product recall.

How do I review the lot number and read the format?

- Look at the surgical gown packaging example picture below:
 - Material number, pointed out by black arrow, is located in the upper right corner (9515)
 - Lot number, pointed out by red arrow, is located in the bottom right corner (17MWQ001)
- In Table 1, you will see the corresponding affected lot number formats for each material number. Please note that the structure of the lot number varies in total characters for different SKUs.
- **Affected lot number format:** X = any character. The characters (X's) will vary, however, for the affected lots, the placement of the letter J will appear in the same position within the sequence of characters as identified in Table 1.
 - For the 9515 example, the affected lot number format for 9515 = XXXXJXXX (i.e. 1234J123, 0000J000, etc.). As long as there are 4 characters before J and 3 characters after, the lot is affected.
 - In the example picture below, the lot number (17MWQ001) does not reflect an affected lot, due to the absence of “J” in relative position of 9515’s lot numbering sequence. Therefore, the customer could continue using this product.



Questions/feedback

GMB-Presource-Notifications@cardinalhealth.com

TABLE 1		
Material	Description	Affected Lot Number Format
90270NB	AAMI 3 NONREIN SCRBNURGOWN BNS L	XXXXXJ or XXXXXJX
90280NC	AAMI 3 ROYALSILK SCRBNURGOWN BNS L	XXXXXJ or XXXXXJX
90370NB	AAMI 3 NONREIN SCRBNURGOWN BNS XL	XXXXXJ or XXXXXJX
90380NC	AAMI 3 ROYALSILK SCRBNURGOWN BNS XL	XXXXXJ or XXXXXJX
90970NA	AAMI 3 FABREINF SCRUB SURG GWN L NS	XXXXXJ or XXXXXJX
91070NB	AAMI 3 FABRNF SCRUB SURG GWN XL NS	XXXXXJ or XXXXXJX
9511	AAMI 3 FABREINF SURG GOWN LARGE	XXXXJXXX
9511NB	AAMI 3 FABREINF SURG GOWN L NS	XXXXXJ or XXXXXJX
9515	GOWN SURG LG W/TOWEL ST BACK	XXXXJXXX
9515NB	AAMI 3 Non-reinf Surg Gown L BNS	XXXXXJ or XXXXXJX
9518	ROYALSILK LARGE	XXXXJXXX
9518NC	AAMI 3 ROYALSILK SURGGOWN BNS L	XXXXXJ or XXXXXJX
9541	AAMI 3 FABRIC REINF SURG GOWN XL	XXXXJXXX
9541NA	AAMI 3 FABRIC REINF SURG GOWN XL NS	XXXXXJ or XXXXXJX
9545	AAMI 3 NON-REINF SURGICAL GOWN XL	XXXXJXXX
9545NC	AAMI 3 NONREIN SURGGOWN BNS XL	XXXXXJ or XXXXXJX
9548	ROYALSILK XLARGE	XXXXJXXX
9548NC	AAMI 3 ROYALSILK SURGGOWN BNS XL	XXXXXJ or XXXXXJX
9571	AAMI 3 FABRIC REINF SURG GOWN 2XL	XXXXJXXX
9571NA	AAMI 3 FABRIC REINFORCED 2XL NS	XXXXXJ or XXXXXJX
9575	AAMI 3 NON-REINF SURGICAL GOWN XXL	XXXXJXXX
9575EL	AAMI 3 NONREINF SURG GOWN 2XL XLONG	XXXXJXXX
9575ELNA	AAMI 3 NONRNF SURG GWN 2XL XLNG BNS	XXXXXJ or XXXXXJX
9575NB	AAMI 3 NON-REINF SURG GOWN 2XL BNS	XXXXXJ or XXXXXJX
9578	ROYALSILK 2XLARGE	XXXXJXXX
9578NB	AAMI 3 ROYALSILK SURGGOWN BNS 2XL	XXXXXJ or XXXXXJX
9585	AAMI 3 NONREINF SURGICAL GOWN 3XL	XXXXJXXX
9585N	AAMI 3 NONREINF SURG GOWN 3XL BNS	XXXXXJ or XXXXXJX
9599	AAMI 3 NONREINF SURGICAL GOWN 4XL	XXXXJXXX
95995	AAMI 3 NON-REINF SURG GWN 3XL XLONG	XXXXJXXX
95995N	AAMI 3 NON-RNF SURG GOWN 3XL LONG	XXXXXJ or XXXXXJX
95998	ROYALSILK 3XLARGE/XLONG	XXXXJXXX
95998N	ROYALSILK 3XLARGE/XLONG NS	XXXXXJ or XXXXXJX
ASG9511	AAMI 3 FABRNF SURG GOWN L 2 TOWELS	XXXXJXXX
ASG9515	AAMI 3 NON-REINF SURG GOWN LG 2 TWL	XXXXJXXX
ASG9541	AAMI 3 FABRNF SURG GWN XL 2 TOWELS	XXXXJXXX
ASG9541EL	AAMI 3 FABRNF SURGGWN XL LONG 2 TWL	XXXXJXXX
ASG9545	AAMI 3 NON-REINF SURG GWN XL 2 TWL	XXXXJXXX