

TEST COMPARISON

	BinaxNOW	QuickVue	iHealth	FlowFlex
CAN IT BE USED AT HOME?	Yes	Yes	Yes	Yes
STORAGE TEMP	36 - 86 F	59 - 86 F	36 - 86 F	36 - 86 F
RESULT TIME	~15 min	~10 min	~15 min	~15 min
WHEN IS PCR CONFIRMATION NEEDED?	Symptomatic NEG & Asymptomatic POS	Symptomatic NEG & Asymptomatic POS	Symptomatic NEG & Asymptomatic POS	Symptomatic NEG & Asymptomatic POS

NOTE: OTC tests are *not* clinical supplies. Tests must be administered by the instruction insert, and patients must be the ones to read results.